


# Safety in numbers



**1** The Blade Remover & Blade Removal Unit have been designed to function with all standard No.3 & No.4 fitment scalpel blades and handles.

Combining simplicity of operation with safety during use these devices assist those who are less confident blade users and especially those working in areas where a sharps disposal facility is not immediately available.

**2** Swann-Morton Disposable Scalpels with retractable blade can assist in reducing cross infection and exposure to blood-borne pathogens arising from accidental sharps injuries and supports existing safe passing protocols.

This simple innovation represents a major breakthrough in handling safety with no compromise in blade quality or control, plus making disposal simpler and safer.

**Swann-Morton®**

**Made in Sheffield - Precisely**

Owlerton Green, Sheffield S6 2BJ

Telephone: 0114 2344231, Sales 0114 2344223

[uksales@swann-morton.com](mailto:uksales@swann-morton.com), [exportsales@swann-morton.com](mailto:exportsales@swann-morton.com)

[www.swann-morton.com](http://www.swann-morton.com)