

Swann-Morton®

SURGICAL KNIFE HANDLES

All handles are individually wrapped. Surgical Handle No's 3, 3L, 5, 7, 9 and B3 take Standard Surgical Blade No's 6, 9, 10, 10A, 11, 11P, E/11, 12, 12B, 12D, 14, 15, 15A, 15C, 15T, D/15 and 16.

Surgical Handle No's 4, 4L and 6 take Standard Surgical Blade No's 18, 19, 20, 21, 22, 22A, 23, 24, 25, 25A, 26, 27 and 36.

The stainless steel B3 Handle, easily gripped, is light, adaptable, autoclavable, and supplied in individual plastic wallets.

With the exception of No's 5 and 6, all handles are also available in Stainless Steel, as well as in nickel. Handle B3 is only available in Stainless Steel.

SMALL FITMENT

LARGE FITMENT

CLEANING



Care must be taken when cleaning handles to avoid damage to the handle fitment, which could result in the blade not fitting correctly.



FOR CLEANING, AND FITTING INSTRUCTIONS - PLEASE SEE OVERLEAF