




Swann-Morton®

Made in Sheffield - *Precisely*


THE WORLD'S LEADING SURGICAL BLADE MANUFACTURER

Surgical Blades for Chiropody and Podiatry. Swann-Morton manufacture a range of carbon or stainless steel sterile surgical blades ideally suited for chiropodists and podiatrists.

Standard Surgical Blades and Handles - No 3 fitment



'Fine' Blades and Handles



PD Blades and Handles



Standard Surgical Blades and Handles - No 3 fitment

Surgical Blade No's 6, 9,10, 10A, 11, 11P, E11, Sabre E11, 12, 12D, 13, 14, 15, 15A, 15C, 15T, Sabre D/15, 16, 40 fit No's 3, 3L, 5, 7, B3 & B3L Handles.

Handles are available in nickel, with selected sizes available in stainless steel. All handles are individually wrapped. The easily gripped B3 is a light, adaptable stainless handle, supplied individually wrapped.

'Fine' Blades and Handles

The 'Fine' range has small, slim, precise blades in stainless steel for delicate operations. Six special handles, including the SF13 and SF23 with independently screwed top sections to facilitate blade-changing, take all blades, supplied in individual metal foil packets, available in dispenser units of 25. Handles are supplied individually in plastic wallets.

PD Blades and Handles

Non sterile, stainless steel blades individually sealed in a view pack and boxed in units of 5 pieces. Stainless steel sterile blades, individually sealed in non autoclavable foil, boxed in units of 25 pieces.

Slim barrel handles are available in stainless steel, with a length of 108mm. An independently screwed top section makes blade changing easier.



All medical devices are CE marked in accordance with the Medical Device Directive (93/42/EEC)



ISO 13485

Blades in carbon steel (sterile/non-sterile) or stainless steel (sterile). Sterile blades individually packed in foil non-sterile individually wrapped in protective paper in foil packets of 5 or in single-peel packets, all boxed in 100's.

Swann-Morton®

Made in Sheffield - *Precisely*

THE WORLD'S LEADING SURGICAL BLADE MANUFACTURER

Enquiry form

How to order

- 1 Enter the type and quantities required in the relevant boxes
- 2 Complete your contact details
- 3 FaxBack the completed form to your usual supplier.

	Type and Quantity	Type and Quantity	Type and Quantity	Type and Quantity	Type and Quantity	Type and Quantity	Type and Quantity
Handle							
Sterile Blade							
Non-Sterile Blade							
Sterile Blade							
Non-Sterile Blade							

Name _____

Company _____

Address _____

_____ Post code _____

Tel. _____ Fax. _____

Email. _____

If you receive this in error, please notify the sender



All medical devices are CE marked in accordance with the Medical Device Directive (93/42/EEC)



ISO 13485

Swann-Morton®

Made in Sheffield - *Precisely*

Owlerton Green, Sheffield S6 2BJ

Telephone: 0114 2344231, Sales 0114 2344223

uksales@swann-morton.com, exportsales@swann-morton.com

www.swann-morton.com