


Small and yet perfectly formed




Silver's, designed specifically for taking small and intricate grafts.

"New" to the range the Silver's stainless steel "miniature" skin graft knife was originally developed by H.L.Silver in Toronto, Canada.

Combined with our surgical stainless double edged razor blades it is ideal for taking small intricate grafts from areas which would not be accessible with the larger standard knives.

Although used mainly within Plastic, Reconstructive and Oral surgery the Silver's can also be used for the debridement of wounds.

The Silver's Skin Graft Knives are individually boxed whilst the blades are sealed in foil packets, sterilised and boxed in units of 10 pieces.



Swann-Morton[®]
Sheffield - England

Swann-Morton Ltd. Penn Works, Owlerton Green, Sheffield S6 2BJ
Telephone: 0114 2344231, Sales 0114 2344223, Fax: 0114 2314966
uksales@swann-morton.com, exportsales@swann-morton.com
orders@swann-morton.com, www.swann-morton.com