



Swann-Morton®

Made in Sheffield - *Precisely*

THE WORLD'S LEADING SURGICAL BLADE MANUFACTURER

Surgical Blades used in Veterinary Surgery. Swann-Morton manufacture a range of carbon or stainless steel sterile surgical blades ideally suited for veterinary surgery.

Standard Surgical Blades and Handles - No 3 fitment



Standard Surgical Blades and Handles - No 4 fitment



Disposable Scalpels



Swann-Morton surgical blades can be used in the surgery for the following:

Small Animals

• Orthopaedics	No.15, No.20	• Hip and Elbow Dysplasia	No.15, No.20
• General Endoscopy	No.15	• Cranial Cruciate Ligament Injuries	No.11, No.11P, No.15, No.20
• Soft Tissue	No.15, 15C, No.20	• Meniscectomy	No.11, No.11P
• Neurosurgery	No.15, No.20	• Bone Fractures and Muscle Diseases	No.15, No.20

Large Animals

• Skin Incisions	No.22	• Animal Pathology PM40 Post Morton Range	• Vaginal Surgery (Horses) 3L, 4L Handles
• Soft Tissue	No.20	• Recto Vaginal Hernias (Horses) 3, 4L Handles	

Specialist Surgical Applications

Standard Surgical Blades and Handles - No 3 fitment

Surgical Blade No's 6, 9,10, 10A, 11, 11P, E11, Sabre E11, 12, 12D, 13, 14, 15, 15A, 15C, 15T, Sabre D/15, 16, 40 fit No's 3, 3L, 5, 7, B3 & B3L Handles.

Handles are available in nickel, with selected sizes available in stainless steel. All handles are individually wrapped. The easily gripped B3 is a light, adaptable stainless handle, supplied individually wrapped.

Standard Surgical Blades and Handles - No 4 fitment

Surgical Blade No's 18, 19, 20, 21, 22, 22A, 23, 24, 25, 25A, 26, 27, 36 fit No's 4, 4L and 6 handles.


Handles are available in nickel, with selected sizes available in stainless steel. All handles are individually wrapped.

Disposable Scalpels

Available with all standard stainless steel blade shapes, complete with plastic handles, individually wrapped in sterile packs of 10. Also available non-sterile with polythene blade guard.

Post Mortem Knife - PM40 and PM40B

Post Mortem Knife with a choice of heavy-duty carbon steel blades. Screw-type fitting to stainless steel handle. Spare blades, individually wrapped, boxed in 10's.



Blades in carbon steel (sterile/non-sterile) or stainless steel (sterile). Sterile blades individually packed in foil non-sterile individually wrapped in protective paper in foil packets of 5 or in single-peel packets, all boxed in 100's.

Swann-Morton®

Made in Sheffield - *Precisely*

THE WORLD'S LEADING SURGICAL BLADE MANUFACTURER

Enquiry form

How to order

- 1 Enter the type and quantities required in the relevant boxes
- 2 Complete your contact details
- 3 FaxBack the completed form to your usual supplier.

	Type and Quantity	Type and Quantity	Type and Quantity	Type and Quantity	Type and Quantity	Type and Quantity	Type and Quantity
Handle							
Sterile Blade							
Non-Sterile Blade							
Sterile Blade							
Non-Sterile Blade							

Name _____

Company _____

Address _____

_____ Post code _____

Tel. _____ Fax. _____

Email. _____

If you receive this in error, please notify the sender



All medical devices are CE marked in accordance with the Medical Device Directive (93/42/EEC)



ISO 13485

Swann-Morton®

Made in Sheffield - *Precisely*

Owlerton Green, Sheffield S6 2BJ

Telephone: 0114 2344231, Sales 0114 2344223

uksales@swann-morton.com, exportsales@swann-morton.com

www.swann-morton.com