

# More comfort less fatigue



## The Swann-Morton range of Ergonomic handles - designed for comfort of grip whilst minimising fatigue during use

The “NEW” range of Swann-Morton Ergonomic handles are designed to maximise performance and feature a textured surface finish for ease and comfort of grip whilst minimising fatigue.

These ergonomically designed handles are used in a Histology, Pathology to Autopsy and Laboratory settings and can be autoclaved to avoid any cross contamination of samples. The handles are also colour coded for ease of recognition.



# Swann-Morton®

ONE OF THE WORLD'S LEADING SURGICAL BLADE MANUFACTURERS

The Black No.5 Ergonomic Handle (1003) is compatible with any of the blade shapes included in the No.3 fitment range such as the 10 and 15 whilst the Maroon No.6 Ergonomic Handle (1004) is compatible with any of the blade shapes included in the No.4 fitment range such as the 22 and 24. Finally there is the Grey No.8 Ergonomic Handle (2618) which is compatible with the pathology/autopsy blades PM60 and PM60 bull nose.



No.5 Ergonomic Handle

No.6 Ergonomic Handle

No.8 Ergonomic Handle

Product	Colour	Product No	Unit	Product	Colour	Product No	Unit
<b>Ergonomic Handles - Non Sterile</b>							
<b>Ergonomic Handle</b>				<b>Ergonomic Handle</b>			
Polypropolyne / Stainless	Black No.5	1003	Each	Polypropolyne / Stainless	Grey No.8	2618	Each
Polypropolyne / Stainless	Maroon No.6	1004	Each				

Swann-Morton and the 'ring pattern' logo are registered trademarks of Swann-Morton Limited and related companies. The products shown are manufactured for use by trained health care professionals and can cause injury if used or handled incorrectly. The information and illustrations used in connection with the products shown on this leaflet are for general guidance only and users must seek proper training and instruction in their use before using or handling such products. Swann-Morton Limited exclude all liability for injury to the extent permitted by law.