



















Symbol	Description / Meaning
	Medical device (Indicates the item is a medical device)
	Sterilized using irradiation (Indicates a medical device that has been sterilized using irradiation)
	Sterilized using ethylene oxide (Indicates a medical device that has been sterilized using ethylene oxide)
	Single Sterile Barrier System (Indicates a single sterile barrier system. This symbol shall be placed adjacent or in combination with the Sterile R, Sterile EO symbols as identified above)
 EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague The Netherlands	Authorized representative in the European Community/European Union (Indicates the authorized representative in the European Community/European Union)
 MEDENVOY SWITZERLAND Gotthardstrassa 28 6302 Zug Switzerland	Authorized representative in Switzerland (Indicates the authorized representative in Switzerland)
	Do not re-use (Indicates a medical device that is intended for one single use)
	Catalogue number (Indicates the manufacturer's catalogue number so that the medical device can be identified)
	Batch code (Indicates the manufacturer's batch code so that the batch or lot can be identified)
	Use-by date (Indicates the date after which the medical device is not to be used)
	Do not use if pack is damaged (Indicates a medical device that should not be used if the package has been damaged or opened)
	Manufacturer (Is accompanied by the name and address of the manufacturer)
	Date of manufacture (Indicates the date when the medical device was manufactured)
	Consult instructions for use or consult electronic instructions for use (Indicates the need for the user to consult the instructions for use)
	Non-sterile (Indicates a medical device that has not been subjected to a sterilization process)
	Caution (To indicate caution is necessary when operating the device, or needs operator awareness in order to avoid undesirable consequences)
	Lock
	Unlock