


















Symbol	Description / Meaning
	Contains hazardous substances (Indicates a medical device that contains substances that can be carcinogenic, mutagenic or reprotox (CMR))
	Medical device (Indicates the item is a medical device)
	Sterilized using irradiation (Indicates a medical device that has been sterilized using irradiation)
	Sterilized using ethylene oxide (Indicates a medical device that has been sterilized using ethylene oxide)
	Authorized representative in the European Community/European Union (Indicates the authorized representative in the European Community/European Union)
	Do not re-use (Indicates a medical device that is intended for one single use)
	Catalogue number (Indicates the manufacturer's catalogue number so that the medical device can be identified)
	Batch code (Indicates the manufacturer's batch code so that the batch or lot can be identified)
	Use-by date (Indicates the date after which the medical device is not to be used)
	Do not use if pack is damaged (Indicates a medical device that should not be used if the package has been damaged or opened)
	Manufacturer (Is accompanied by the name and address of the manufacturer)
	Date of manufacture (Indicates the date when the medical device was manufactured)
	Consult instructions for use or consult electronic instructions for use (Indicates the need for the user to consult the instructions for use)
	Non-sterile (Indicates a medical device that has not been subjected to a sterilization process)
	Caution (To indicate caution is necessary when operating the device, or needs operator awareness in order to avoid undesirable consequences)
	Lock
	Unlock