




# Finger on the Trigger

## An extended Retractable Scalpel range with Sharps Safety compliance in mind.


The EU Council Directive 2010/32/EU recommends that medical devices incorporating engineered safety mechanisms and protection should be “considered” by staff as “part” of their risk assessment process. This would include devices such as the Swann-Morton Retractable Scalpel now extended to include a Safety Stitch Cutter.



**1** Temporary lock in the forward position has been strengthened so that the blade will stay extended even if it encounters a hard surface or bone during use.



**2** The slider mechanism now contains a notch in the retracted position that holds the blade more securely in the passing/retracted position.



**3** A permanent lock position has been added to the retracted position. This requires a conscious amount of force to engage the slider into the permanent lock position.




# Swann-Morton®

THE WORLD'S LEADING SURGICAL BLADE MANUFACTURER

## Retractable Disposable Scalpel

Designed for single patient use only, disposable, sterile or non-sterile and Latex safe.



### Retractable Disposable Scalpel - Sterile

Plastic Handle fitted with Stainless Steel Surgical Blade, individually packed

Product	Product No	Unit
<b>25 scalpels individually packed</b>		
10	3901	25
10A	N/A	
11P	3991	25
15	3905	25
20	3906	25
21	3907	25
22	3908	25
23	3910	25
24	3911	25
25A	3915	25
Stitch Cutter	3926	25


Other Blade shapes available on request

### Retractable Disposable Scalpel - Non-Sterile

Plastic Handle fitted with Stainless Steel Surgical Blade, bulk packed

Product	Product No	Unit
<b>500 scalpels bulk packed</b>		
10	4001	Box of 500
10A	4002	Box of 500
11P	4091	Box of 500
15	4005	Box of 500
20	4006	Box of 500
21	4007	Box of 500
22	4008	Box of 500
23	4010	Box of 500
24	4011	Box of 500
25A	4015	Box of 500
Stitch Cutter	4026	Box of 500

Other Blade shapes available on request



Swann-Morton and the 'ring pattern' logo are registered trademarks of Swann-Morton Limited and related companies. The products shown are manufactured for use by trained health care professionals and can cause injury if used or handled incorrectly. The information and illustrations used in connection with the products shown on this leaflet are for general guidance only and users must seek proper training and instruction in their use before using or handling such products. Swann-Morton Limited exclude all liability for injury to the extent permitted by law.

**Swann-Morton®**

Sheffield - England

Swann-Morton Ltd. Penn Works, Owlerton Green, Sheffield S6 2BJ  
 Telephone: 0114 2344231, Sales 0114 2344223, Fax: 0114 2314966  
 uksales@swann-morton.com, exportsales@swann-morton.com  
 orders@swann-morton.com, www.swann-morton.com

RET/issue 3 12/14