

# The Pathologists' first choice



## A versatile and comprehensive range for Pathologists

As part of their diverse range of Surgical Blades and Handles, Swann-Morton offer two handle and blade combinations specifically for use within Post Mortem examination.

The PM40 and PM60 range are used by Pathologists during Forensic and Clinical Autopsy procedures and for Post Mortem diagnosis, the quality of which meets the same exacting high standards for which Swann-Morton are renowned amongst the medical and related health care professions worldwide.

The PM40B and PM60B are "bull nose" blades for added safety and protection against accidental sharps injuries.



**Swann-Morton®**  
Made in Sheffield - *Precisely*




FM 73368

# Swann-Morton®

Made in Sheffield - *Precisely*

THE WORLD'S LEADING SURGICAL BLADE MANUFACTURER

The heavy duty carbon steel blade No's PM40 and PM40B are both used during Post Mortem examinations and fit the standard "blue" coated or "plated" PM40 Handles using a traditional screw fitting.



An autoclavable, solid-cast stainless steel handle was more recently added to the range with a screw fitting incorporated in the top section only. This handle is supplied singly in a protective plastic wallet. The length of the handle with blade is approximately 220mm.

Alternatively we have the PM8 handle which is a light yet strong nickel handle with the standard surgical range type of fitment\*. This is again supplied singly in a protective plastic wallet for use with the carbon steel blade No's PM60 and PM60B\*\*. The length of the handle with blade is approximately 185mm and the PM60 range is the preferred choice of Pathologists and Morticians in the USA.



PM40 Range blades			PM60 Range blades		
Product	Product No	Unit	Product	Product No	Unit
PM40 Standard	2551	10 blades	PM60 Standard	2560	50 blades
PM40 B	2552	individually packed	PM60 B	2561	individually packed
PM40 Range Handles			PM60 Range Handles		
Product	Product No	Unit	Product	Product No	Unit
PM40 Stainless	2650	1 handle	PM8	2608	1 handle
PM40 Standard	2652	1 handle	For use with 60 and 60 B blades		
PM40 Plated	2653	1 handle			

\* The blade fitment on the PM8 Handle is "not" compatible with other standard range surgical blades. \*\* The PM60B blade can also be referred to as the PM70 blade. Note: A number of blades from the heavy duty Swann-Major blade and handle range plus the No.22 standard blade shape are utilised by both Pathologists and Morticians. For further related product information please visit [www.swann-morton.com](http://www.swann-morton.com)

Swann-Morton and the 'ring pattern' logo are registered trademarks of Swann-Morton Limited and related companies. The products shown are manufactured for use by trained health care professionals and can cause injury if used or handled incorrectly. The information and illustrations used in connection with the products shown on this leaflet are for general guidance only and users must seek proper training and instruction in their use before using or handling such products. Swann-Morton Limited exclude all liability for injury to the extent permitted by law.